



4. How to launch DeepKymoTracker

(User guide)


Go to command line and start typing **Anaconda Navigator**. In the open menu, click on the icon **Anaconda Navigator**.




On the top of the screen, you will see a small window with words base(root). Click on the arrow on the right, and then click on **py36** in the dropdown menu.




Here is what you get. All that is left for you to do is to click **Launch** button under Spyder icon.



This is the Spyder window. Open file **GUI_launch.py** (it is inside DeepKymoTracker folder; click on the second left icon to open file menu and navigate to this file) and launch it by clicking **green arrow icon** on the top of the screen.



If nothing happens, press the **button with the icon** at the bottom of the screen as shown in the picture below:



Now, you can see the title page of DeepKymoTracker. It has 5 more pages (windows) each devoted to a certain step of the pipeline. You can navigate to the page of interest either by using button NEXT at the bottom or by pushing one of the green buttons GO TO STEP -1, GO TO STEP -2, etc.

